

Reversible Tube Racks, 96-WELL Heathrow Scientific

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Research & Development](#) / [Refrigeration and freezing](#) / [Tube boxes & freezer racks](#) / Reversible Tube Racks, 96-WELL Heathrow Scientific

Product Data

Easy-to-grasp tabs for sample handling and storage

96 wells for 0.5 mL tubes on one side and 96 wells for 1.5/2.0 mL tubes on the reverse side

Optimal Function & Identification

- Wells arranged in standard 8 x 12 array
- Molded in grid lines separate 0.5 mL wells making them easier to see
- Raised rings around the 1.5 mL/2.0 mL wells make tubes easier to remove
- Alphanumeric identification for quick location and retrieval of tubes

Secure Lid

- Translucent lid snaps onto rack tabs for secure, visible storage
- Spacious enough to fit most screw-top micro-centrifuge tube
- Stackable, saving valuable space

Specifications

- Wells: 96
- Made of polypropylene
- Tube Sizes: 0.5 mL and 1.5/2.0 mL
- Dimensions: 9.6 x 4.7 x 2 in (24.4 x 12 x 5.1 cm)
- Comes in packs of 5

Item No.	colour	Price
15.2345.48	natural	CHF 93.75
15.2345.49	blue	CHF 93.75
15.2345.50	green	CHF 93.75
15.2345.51	pink	CHF 93.75
15.2345.52	yellow	CHF 93.75
15.2345.53	assorted: blue, green, pink, yellow and orange	CHF 93.75

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 27.04.2026 at 22:54 oclock.

www.huberlab.ch