

Sample Bottle, Wide Mouth Wheaton



Product in eShop

[Home](#) / [Assortment](#) / [Quality Assurance](#) / [Sampling](#) / Sample Bottle, Wide Mouth Wheaton

Product Data

- For use with Sub Surface Grab Sampler® III
- Wide mouth borosilicate glass bottle
- 70-400 screw cap with PTFE liner
- Autoclavable
- These wide mouth sample bottles can be repeatedly autoclaved and come with PTFE-lined caps

Item No.	volume	Price
8.9904.76	1000 ml	CHF 339.95

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 20.04.2026 at 15:59 o'clock.