

# Samplepack Centrifuge Tubes 50 ml Greiner Bio-One

[Product in eShop](#)

[Home](#) / [Brands](#) / [Greiner Bio-One](#) / [Samplepack Centrifuge Tubes 50 ml Greiner Bio-One](#)

## Product Data

The Greiner Bio-One tube range is extremely versatile and meets a wide variety of requirements. The tubes offer excellent thermal, mechanical and chemical stability. They also feature precise graduations and a large white label area.

- Free from detectable DNase, RNase and human DNA
- Free from detectable endotoxins, free from cytotoxic substances
- Polypropylene has high thermal, mechanical and chemical resistance, making it very well suited for storing chemical and biological samples
- With conical bottom for optimal centrifugation
- With graduation and writing field
- With screw cap

**Sample pack contains:** 7.210 261 and 7.227 261-N

Item No.	Price
7.227 999	for free

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 14.05.2026 at 20:51 o'clock.