

## SteriPlast® Kit – sterile sampling set

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Life Science](#) / [Microbiology](#) / [Sampling](#) / SteriPlast® Kit – sterile sampling set

### Product Data

The sampling set SteriPlast® Kit consists of ten sample spatulas or ten sample scoops and ten SteriBag Premium tamper-proof sample bags.

Thus the SteriPlast® Kit is the ideal combination for sampling where purity and sterility is essential.

All parts are sterilised by gamma rays and suitable for use with foodstuffs.

The scoops or spatulas are clean room manufactured and individually packed. The sample bag SteriBag Premium has an integrated sealing and is liquid-tight. A unique tracking code on both the bag and the tear-off field guarantees complete documentation and traceability.

SteriPlast® Kit including scoop 50 ml or spatula 150 mm insertion depth, and SteriBag Premium 300 ml.

- Scoop 50 ml & bag SteriBag Premium 300 ml
- Spatula 150 mm insertion depth & bag SteriBag Premium 300 ml

Item No.	execution	Price
15.5378.92	Scoop & bag	CHF 44.80
15.5378.93	Spatula & bag	CHF 41.95

### More product images



\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 22.05.2026 at 15:30 oclock.

[www.huberlab.ch](http://www.huberlab.ch)