

Test strips, peracetic acid, MQuant™

[Product in eShop](#)[Home](#) / [Brands](#) / [Merck](#) / Test strips, peracetic acid, MQuant™

Product Data

MQuant™ (MQ) test strips can be considered as a mobile laboratory on a few square millimeters of plastic strip. As a screening method they provide the user with a rapid overview of the situation, which helps in making the choice of what, and perhaps if, further analysis is required, thereby leading to potentially substantial savings in time and analytical resources. The strip is a biologically degradable polyester foil and the low reagent content, in the reaction zones, assists problem free waste disposal.

- Easy to handle: just dip, wait and measure visually against colour scale
- Safe: no auxiliary reagents necessary
- Immediate results within 2 minutes enables fast decisions
- Portable: can be used inside or outside the laboratory

Item No.	execution	Price
1.10001.0001	Range: 100-150-200-250-300-400-500 mg/l CHCOH	CHF 40.08 CHF 134.35
1.10084.0001	Range: 5-10-20-30-50 mg/l CHCOH	CHF 91.50

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 19.04.2026 at 13:37 o'clock.