

Bosotherm flex Fieberthermometer

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Medicine & Healthcare](#) / [Personal care](#) / [Thermometers](#) /

Product Data

Das bosotherm flex Fieberthermometer für die klassische Fiebermessung ermittelt die Körpertemperatur in nur 10 Sekunden und verfügt über eine flexible Spitze für eine möglichst sanfte Messung.

Produktdetails:

- Schnelle rektale Fiebermessung in nur 10 sec
- Flexible Spitze für eine sanfte Fiebermessung
- Antiallergische Goldbeschichtung, 24 Karat vergoldet
- 100% wasserdicht
- Speicher für den letzten Messwert
- Praktisches Klarsichtetui
- inklusive Batterie
- 2 Jahre Garantie

Item No.	Price
3.7666.05	CHF 10.05

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 19.05.2026 at 07:15 o'clock.