

## neodisher® IR



Product in eShop

Home / [Assortment](#) / [Safety & Hygiene](#) / [Cleaning](#) / [Cleansers & soaps](#) / neodisher® IR

### Product Data

For the thorough cleaning of surgical instruments made of stainless steel in immersion baths and ultrasonic baths. Tarnish, rust film and extraneous rust are removed. Observe product data sheet before use.

- UN 1805 PHOSPHORSÄURELÖSUNG, 8, III, E

Item No.	execution	capacity	Price
18.9420.41	bottle	12x 1 liter	CHF 302.30
18.9420.42	bottle	1 liter	CHF 25.20

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 27.04.2026 at 22:53 o'clock.