

PLR- Drehzahlfeinreguliertventil Zubehör zu PLR- Druckluft-Laborrührer



Product in eShop

Home / [Assortment](#) / [Research & Development](#) / [Synthesis](#) / [Stirrers](#) /

Product Data

| Item No. | description | Price |
|------------|---|------------|
| 17.1266.01 | PLR- Drehzahlfeinreguliertventil | CHF 209.15 |
| 17.1266.02 | PLR - Fingerdruck-Oeler, für Normale-und Tacho- Version | CHF 257.60 |
| 17.1266.30 | PLR-Plattenstativ BxTxL 210x315x1000 mm | CHF 308.30 |
| 17.1266.31 | PLR-Kreuzklemme 16 mm, drehbar | CHF 96.60 |
| 17.1266.33 | PLR-Verbindungskupplung VK 10 x 8 | CHF 119.00 |

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 23.04.2026 at 15:18 o'clock.