

Compressor-cooled incubator ICP / ICPeco Memmert

[Product in eShop](#)

Home / [HUBERLAB.green](#) / Compressor-cooled incubator ICP / ICPeco Memmert

Product Data

Incubate, dry and store in the cooled incubator ICPeco.

This environmentally-friendly incubator is cooled with climate-friendly CO₂. Thanks to the excellent thermodynamic properties of the refrigerant CO₂ (R744) and the finely adjusted control technology, the ICPeco cooled incubator is both powerful and high-precision. Without critical temperature overshoots, it keeps the temperatures exactly at the setpoint.

- temperature range up to +60 °C
- 4 model sizes (108 to 749 litres volume)
- double doors standard for all models: Prevention of contamination and drops in temperature, and at the same time an optimum view of the sensitive load through wide-area interior glass doors
- very good air circulation in the working chamber through sophisticated ventilation technology
- intelligent defrosting function

Item No.	model	vol. l	Price
17.1180.05	ICP110	108 l	CHF 8'992.50
17.1180.06	ICP260	256 l	CHF 10'278.00
17.1180.07	ICP450	449 l	CHF 14'433.75
17.1180.08	ICP750	749 l	CHF 16'953.75
17.1180.10	ICP260eco	256 l	CHF 11'165.00
17.1180.11	ICP450eco	449 l	CHF 15'306.00

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 21.04.2026 at 06:37 o'clock.

www.huberlab.ch