

Pioneer™ PX Analytical and Precision Balances Ohaus



Product in eShop

[Home](#) / [Assortment](#) / [Research & Development](#) / [Weighing](#) / [Scales](#) / Pioneer™ PX Analytical and Precision Balances Ohaus

Product Data

High-performance, economical balances with high accuracy and repeatability for applications in laboratories, industry and education. For basic weighing, parts counting, percent weighing, dynamic weighing and density determination.

- 2 line, backlit LCD display
- 4 quick keys enable the most common weighing modes
- USB and RS 232 interface for data exchange with PCs and PLCs
- User-selectable span calibration points
- Real time clock with GLP/GMP data
- Cast metal lower housing, removable stainless steel weighing pan
- Integrated weigh-below-hook
- Auto-tare
- Many weighing units selectable, a custom unit can be defined
- Overload, underload and stability indicator
- 4 menu languages

Please note:

Models with extra code /E with external calibration, all other models with internal calibration, models with extra code M are approved instruments according to EN 45501, OIML R76.

Specifications

Permissible ambient temperature: 10 ... 30 °C

Weighing plate (Ø): 180 mm

Height in draftshield: 235 mm

Dimensions (W x D x H): 209 x 321 x 98 mm

Weight: 3.5 kg

Power supply: 100 ... 240 V, 50/60 Hz

Item No.	type no.	weighing range	display precision	pan	Price
----------	----------	----------------	-------------------	-----	-------

16.4298.03	PX124 analytical balance	120 g	0.0001 g	90 mm	CHF 1'765.00
16.4298.04	PX224 analytical balance	220 g	0.0001 g	90 mm	CHF 1'900.00
16.4298.05	PX323 precision balance	320 g	0.001 g	120 mm	CHF 1'345.00
16.4298.09	PX3202/E precision balance	3200 g	0.01 g	180 mm	CHF 1'215.00
16.4298.10	PX3202M precision balance	3200 g	0.01 g	180 mm	CHF 1'775.00
16.4298.11	PX224/E analytical balance	220 g	0.0001 g	90 mm	CHF 1'685.00
16.4298.12	PX4201 precision balance	4200 g	0.1 g	180 mm	CHF 1'540.00
16.4298.13	PX4201/E precision balance	4200 g	0.1 g	180 mm	CHF 1'145.00
16.4298.14	PX4201M precision balance	4200 g	0.1 g	180 mm	CHF 1'675.00
16.4298.15	PX6202 precision balance	6200 g	0.01 g	180 mm	CHF 1'645.00
16.4298.16	PX6202/E precision balance	6200 g	0.01 g	180 mm	CHF 1'380.00
16.4298.17	X6202M precision balance	6200 g	0.01 g	180 mm	CHF 1'925.00
16.4781.97	PX225D Semi-Micro	82 g/ 220 g	0.00001/0.0001 g	80 mm	CHF 4'550.00
16.4853.40	PX85 Semi-Micro	82 g	0.00001 g	80 mm	CHF 3'590.00

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 09.05.2026 at 04:43 o'clock.