

## Titrating burette automatic Bürkle

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Research & Development](#) / [Titration](#) / Titrating burette automatic Bürkle

### Product Data

Burette with Schellbach stripes type Schilling / Bürkle with original Bürkle dosing armature. A universal aid in the laboratory, electroplating industry, water treatment plants, for field analyses, etc.

Precise dosing using the press button, exact fine titration with the micro-screw and precise discharge jet. Graduation according to Category B (Tolerances according to DIN EN ISO 385).

- Adjusted for Ex
- Burette made of borosilicate glass 3.3
- Resistant to chemicals
- Blue ring graduations (Colordur)
- Left or right hand operation possible
- Sturdy stand
- Precise discharge jet, allowing very fine drops
- Parts in contact with the medium: PE, PP, natural rubber, borosilicate glass

Item No.	volume	tol. +/-	storage bottle	division	Price
1.9695.05	5 ml	0,015 ml	500 ml	0.05 ml	CHF 75.55
1.9695.10	10 ml	0,030 ml	500 ml	0.05 ml	CHF 75.55
1.9695.15	15 ml	0,045 ml	500 ml	0.10 ml	CHF 75.55
1.9695.25	25 ml	0,045 ml	1000 ml	0.10 ml	CHF 75.55
1.9695.50	50 ml	0,045 ml	1000 ml	0.10 ml	CHF 88.00

[More product images](#)



\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 17.04.2026 at 08:08 o'clock.

www.huberlab.ch