


Diagnostic slides Diagxer 385 W

[Product in eShop](#)


[Home](#) / [Assortment](#) / [Medicine & Healthcare](#) / [Microscopy](#) / [Slides](#) / Diagnostic slides Diagxer 385 W

Product Data

Dimensions 75 x 25 mm, ground edges 90°

Item No.	execution	colour	number of reaction fields	d of reaction field	Price
10.0230.53	Epoxy coated	black	10	6,7 mm	CHF 121.60
10.0231.11	Teflon coating	white	1 (oval)	---	CHF 220.20
10.0232.03	PTFE-coating	black	3	14 mm	CHF 266.00
10.0232.10	PTFE-coating	black	10	6,7 mm	CHF 333.75

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 02.05.2026 at 17:53 o'clock.