

Slide boxes antistatic



Product in eShop

[Home](#) / [Assortment](#) / [Medicine & Healthcare](#) / [Microscopy](#) / [Slide magazines](#) / [Slide boxes antistatic](#)

Product Data

- Black
- Numbered plastic rack, antistatic, resistant to xylene
- With numbered index in the cover
- For micro slides 76 x 26 mm

With wooden frame:

- Bottom part and lid are connected by metal hinges and clasps

Item No.	for number of slides	description	Price
10.2702.25	25	with wooden frame	CHF 20.50
10.2702.50	50	with wooden frame	CHF 27.40
10.2702.99	100	with frame of PP	CHF 27.70

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 14.05.2026 at 20:54 oclock.

www.huberlab.ch