


Tape thermometer testoterm®, irreversible, 8 temperatures testo

[Product in eShop](#)

Home / [Labortops Q2 2026](#) / Tape thermometer testoterm®, irreversible, 8 temperatures testo

Product Data

Strips for 8 temperatures. 5 models cover the measuring range from 37 to 260 °C, with each strip showing 8 different temperatures.

Item No.	measuring range °C	Price
4.2406.02	37 - 65 °C	CHF 28.75 CHF 40.00
4.2406.03	71 - 110 °C	CHF 28.75 CHF 40.00
4.2406.04	116 - 154 °C	CHF 28.75 CHF 40.00
4.2406.05	161 - 204 °C	CHF 28.75 CHF 40.00
4.2406.06	204 - 260 °C	CHF 28.75 CHF 40.00

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 04.05.2026 at 03:02 o'clock.