

# Suspension Culture Flasks CELLSTAR® Greiner Bio-One

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Life Science](#) / [Cell Biology](#) / [Cell culture consumables](#) / [Cell culture flasks](#) / [Suspension Culture Flasks CELLSTAR® Greiner Bio-One](#)

## Product Data

Hydrophobic surface, ideal for suspension cultures, hybridoma and embryonic stem cells

- Suspension culture flasks with standard screw cap (without filter)
- White standard screw caps with ventilation position
- Canted neck
- Graduation on both sides
- Sterile and user-friendly packaging

Item No.	vol.	TC treated	packaging	Price
7.658 190	250 ml	-/Suspension	24 packs x 5 pieces	CHF 271.78
7.660 190	550 ml	-/Suspension	10 packs x 5 pieces	CHF 221.55
7.661 190	650 ml	-/Suspension	10 packs x 4 pieces	CHF 177.24
7.690 190	50 ml	-/Suspension	20 packs x 10 pieces	CHF 236.34
7.690 195-TRI	---	---	---	CHF 332.49

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 25.05.2026 at 22:32 o'clock.