

Control seal close-it Bürkle



Product in eShop



[Home](#) / [Assortment](#) / [Laboratory consumables & instruments](#) / [Consumables](#) / [Closure seals & labels](#) / Control seal close-it Bürkle

Product Data

close-itfood is an aluminium backed, vapour-proof adhesive label with extremely strong adhesion. Sacks, boxes etc. which have been pierced with samplers can be hermetically resealed. Can even adhere to slightly moist or dusty sacks coated with fine powders, e.g. plaster, flour etc

- Can be written on with ballpoint pen, pencil or felt-tip marker
- Two sizes: 95 x 95 and 150 x 150 mm

Item No.	format	colour	rolls of	Price
15.5303.01	95 x 95 mm	rot	500 Stk	CHF 79.40
15.5303.02	95 x 95 mm	grün	500 Stk	CHF 79.40
15.5303.03	95 x 95 mm	gelb	500 Stk	CHF 79.40
15.5303.04	95 x 95 mm	blau	500 Stk	CHF 79.40
15.5303.05	95 x 95 mm	schwarz	500 Stk	CHF 79.40
15.5303.06	95 x 95 mm	weiss	500 Stk	CHF 79.40
15.5303.07	150 x 150 mm	rot	250 Stk	CHF 111.25
15.5303.08	150 x 150 mm	grün	250 Stk	CHF 111.25
15.5303.09	150 x 150 mm	gelb	250 Stk	CHF 111.25
15.5303.10	150 x 150 mm	blau	250 Stk	CHF 111.25
15.5303.11	150 x 150 mm	schwarz	250 Stk	CHF 111.25
15.5303.12	150 x 150 mm	weiss	250 Stk	CHF 111.25

Accessories

Item No.	Description	Price
15.5303.13	Close-it dispenser Bürkle	from CHF 124.60

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 24.05.2026 at 11:25 oclock.

www.huberlab.ch