

## Wall Dosing Machines 250 ml

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Laboratory consumables & instruments](#) / [Liquid Handling](#) / [Dose](#) / [Dispenser](#) /

### Product Data

Wand-Dosiergerät für unkompliziertes, sauberes und exaktes Dosieren von Harzen, Lacken, Peroxiden, Epoxidharz sowie vielen anderen pumpfähigen Medien – ohne Verkleben und Verschmutzen. Für den stationären Einsatz am Arbeitsplatz.

- Werkstoffe: Glas, PTFE, PA und Edelstahl
- Einfache Reinigung
- Wandmontage vorbereitet
- Ansaugschlauch PA

Item No.	description	designation	vol.	Price
15.6111.02	Ansaugschaluch 22 mm Ø	kleinste Teilung 5 ml	250 ml	CHF 1'720.32

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 10.05.2026 at 23:26 o'clock.